

Dipolar Couplings in Solids: Measurement and Assignment Using Novel Techniques

**Robert H. Havlin
July 30, 2002**

44th Rocky Mountain Conference

NMR Crystallography ...

The Goal:

Measure enough distances to determine structure in a single experiment.

Liquid State

PrP^C



Many methods including NOE
yield distance constraints for
determining structure

Solid State


PrP^{SC}

Solid


Few methods exist for the
simultaneous determination of
multiple distances

Recoupling Methods


The static experiment contains **everything** including dipolar couplings



MAS improves resolution, however dipolar couplings are **lost**.





Recoupling techniques such as DRAMA, C7,REDOR, RN, etc. **reintroduce** dipolar couplings under MAS.



However, all of these methods for determining dipolar couplings have a **theta** dependence leading to **anisotropic** spectra.

Field Cycling Experiments

Field cycling experiments utilize the fact that in zero field
a dipolar coupling exists as 3 sharp lines:



Zero Field at High Field





$$H_D^{\text{ZF}} = \omega_D \left(\vec{I}^1 \cdot \vec{I}^2 - 3 \left(\vec{I}^1 \cdot \hat{r}_{12} \right) \left(\vec{I}^2 \cdot \hat{r}_{12} \right) \right)$$
$$= \omega_D \sum_{m=-2}^2 (-1)^m A_{2,m}^{12}(\theta_{12}, \phi_{12}) T_{2,-m}^{12}$$

$$H_D^{\text{HF}} = \omega_D A_{2,0}^{12} T_{2,0}^{12}$$

$$= \omega_D \sum_{l=0,2,4} C(2,2,l,0) F_{l,0}$$

$$\overline{H} \propto F_{0,0} = \sigma H_D^{\text{ZF}}$$

Tycko's ZFHF Results



Inspiration

We want an isotropic spectrum similar to zero field.

$$H_D \propto \omega_D (3 \cos^2 \theta - 1)$$

So we have $\cos^2 \theta \dots$

How can θ be eliminated
to leave ω_D ?

Generating an Isotropic Signal

$$\rho(0) = I_z^1 + I_z^2$$

$$H = \omega_D \textcolor{green}{h}(\theta) (3I_x^1 I_x^2 - \overrightarrow{I^1} \overrightarrow{I^2})$$

for a time t not at the magic angle

$$h(\theta) = h \frac{3 \cos^2 \theta - 1}{2}$$

$$\begin{aligned} \rho(t) &= \cos\left(\frac{3}{2}\omega_D \textcolor{green}{h}(\theta)t\right) I_z + i \sin\left(\frac{3}{2}\omega_D \textcolor{green}{h}(\theta)t\right) (T_{2,2} - T_{2,-2}) \\ &= \rho_0(t) + \rho_2(t) \end{aligned}$$

$$\rho_o(t) + \rho_2(t)$$

$$H^{EVO} = \omega_D (\textcolor{blue}{g} A_{2,-2} T_{2,2} + \textcolor{blue}{g} A_{2,2} T_{2,-2})$$

for a time τ

$$g(\theta) = \frac{3}{2} g \sin^2 \theta$$

$$A_{22} = \frac{3}{4} \sin^2 \theta \cdot e^{i2\phi}$$



$$\begin{aligned} \langle I_z(t, \tau) \rangle_1 &= \text{Tr}(\rho_o(t, \tau) I_z) \\ &= \cos(\omega_D \textcolor{green}{h}(\theta) t) \end{aligned}$$

$$\cdot \cos(\omega_D \textcolor{blue}{g}(\theta) \tau)$$

$$\begin{aligned} \langle I_z(t, \tau) \rangle_2 &= \text{Tr}(\rho_2(t, \tau) I_z) \\ &= -\cos(2\phi) \\ &\quad \cdot \sin(\omega_D \textcolor{green}{h}(\theta) t) \\ &\quad \cdot \sin(\omega_D \textcolor{blue}{g}(\theta) \tau) \end{aligned}$$

$$H_{\pm 1}^{DET} = \omega_D \mathbf{f}(A_{2,\pm 1}T_{2,2} - A_{2,\mp 1}T_{2,-2})$$

for a constant time τ_{DET}

$$\begin{aligned} S_1 &= \langle I_z(t, \tau) \rangle_1 = -\cos(\omega_D \mathbf{h}(\theta)t) \\ &\quad \cdot \cos(\omega_D \mathbf{g}(\theta)\tau) \\ &\quad \cdot \sin^2(\omega_D \mathbf{f}(\theta)\tau_{DET}) \end{aligned}$$

$$\begin{aligned} S_2 &= \langle I_z(t, \tau) \rangle_2 = \cos^2(2\phi) \\ &\quad \cdot \sin(\omega_D \mathbf{h}(\theta)t) \\ &\quad \cdot \sin(\omega_D \mathbf{g}(\theta)\tau) \\ &\quad \cdot \sin^2(\omega_D \mathbf{f}(\theta)\tau_{DET}) \end{aligned}$$

using:

$$\cos a \cos b \mp \sin a \sin b = \cos(a \pm b)$$

$$S_1 \pm 2S_2 \propto \cos(\omega_D [\mathbf{h}(\theta)t \pm \mathbf{g}(\theta)\tau])$$

$$\propto \cos \left[\mathbf{h}\omega_D t(3\cos^2(\theta)-1) \pm \mathbf{g}\omega_D \tau \sin^2(\theta) \right] \text{ Final Signal}$$

adjust t and τ to add


$$\propto \cos(\sigma \omega_D t_1)$$

Scaling Factor $\sigma = \frac{3ht}{2(t+\tau)} \omega_D$


$$\sigma_{\max} = 0.30$$

Homonuclear Isotropic Evolution

HOMIE Pulse Sequence




HOMIE Sequence Details




experiment	Hdip phase	DET phase	post-C7	receive phase
S ₁ {	X ₁	0	0	-1
	X ₂	0	90	1
	Y ₁	90	0	-1
	Y ₂	90	90	1
S ₂ {	2x X ₃	0	0	1
	2x X ₄	0	90	-1
	2x Y ₃	90	0	-1
	2x Y ₄	90	90	1

remaining signal $\approx 23\%$

Homonuclear Results




Local Field Spectroscopy



Caldarelli, Hong, Emsley and Pines, JPC 100, 18699.

Linear with the couplings


Heteronuclear Isotropic Sequence




Generate and detect

$$\sin^2 \theta + \cos^2 \theta = 1$$


All while eliminating homonuclear dipolar couplings



$$\overline{H}_C \propto (3 \cos^2 \theta_R - 1)(3 \cos^2 \theta - 1) I_X S_Z$$





$$\overline{H}_S \propto \sin^2 \theta_R \sin^2 \theta (e^{i2\phi} I_+ + e^{-i2\phi} I_-) S_Z$$




$$\overline{H}_{\text{Det}} \propto \sin 2\theta_R \sin 2\theta (e^{i\phi} I_+ - e^{-i\phi} I_-) S_Z$$

Heteronuclear Results

Experiment



Simulation




Scaling factor = 0.08


Future Possibilities

Using selectively introduced labels or even naturally occurring ones, we can approach NMR crystallography.

Additionally, angles can also be determined!



Acknowledgements

A photograph of a large, multi-trunked pine tree in the foreground, its branches reaching across the frame. Behind it is a white, two-story building with a prominent dark steeple or tower. The sky is clear and blue.

**Jamie D. Walls
Wyndham B. Blanton
Prof. Alex Pines
Tanya Mazur
Pines Group**

**NSF and DOE
for Funding**